

Product Datasheet

CD80 antibody (orb607882)

Description

Mouse monoclonal antibody to CD80 (B7-1)

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant fragment of human CD80 protein (around aa 35-142) (exact sequence is proprietary)
Target	CD80
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b, kappa
Clonality	Monoclonal
MW	60kDa
Uniprot ID	P33681
Entrez	941
Dilution Range	ELISA (For coating, order Ab without BSA); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.