

Product Datasheet

CD86 antibody (FITC) (orb607883)

Description

Mouse monoclonal antibody against CD86

	conjugated to FITC
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	ARH-77 (B-lymphoblastoid cell line)
Target	CD86
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab conjugated to FITC
Storage	Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	70kDa
Uniprot ID	P42081
Entrez	942
Dilution Range	Flow Cytometry (1-2ug/million cells); Immunofluorescence (1-2ug/ml); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.