



Product Datasheet

GCN2 (phospho-Thr667) antibody (orb6079)

Description	GCN2 (phospho-Thr667) antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GCN2 around the phosphorylation site of Thr667 (PG(p-T)PP)
Target	phospho-GCN2 (Thr667)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	181 kDa
Uniprot ID	Q9P2K8
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of mouse liver tissue (Dilution at...