

Product Datasheet

HIV1-p24 antibody (CF488) (orb607953)

Description

Mouse monoclonal antibody against Human

Species/Host

Mouse

Conjugation

Unconjugated

Tested

ICC, IF

Applications
Immunogen

Recombinant HIV-1 Gag p24 protein

Preservatives

Antibody Purified from Bioreactor Concentrate by Protein A/G and conjugated to various reporter molecules. Prepared in 10mM PBS with 0.05% rAlbumin and 0.05% azide. Contact us if you require this Ab in a different format.

Concentration

Purified Ab conjugated to CF488

Storage

Antibody with azide - store at 4 to 8°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa

Clonality

Monoclonal

MW

24kDa (mature); 55kDa & 41kDa (precursors)

Uniprot ID

P04585

Entrez

155348

Dilution Range

Immunocytochemistry (Acetone or paraformaldehyde fixed) (1:50-1:100);Immunofluorescence (1:50-1:100);Optimal dilution for a specific application should be determined.

Expiration Date

12 months from date of receipt.