

Product Datasheet

KRT77 & KRT76 antibody (FITC) (orb607964)

Description

Mouse monoclonal antibody against

conjugated to FITC

Species/Host

Mouse

Reactivity

Bovine, Canine, Gallus, Human, Monkey,
Mouse, Rabbit, Rat

Conjugation

Unconjugated

Tested

FC, IF

Applications
Immunogen

Human epidermal keratin

Target

KRT77 & KRT76

Preservatives

Antibody Purified from Bioreactor
Concentrate by Protein A/G and conjugated
to various reporter molecules. Prepared in
10mM PBS with 0.05% rAlbumin and 0.05%
azide. Contact us if you require this Ab in a
different format.

Concentration

Purified Ab conjugated to FITC

Storage

Antibody with azide - store at 4 to 8°C.
Antibody is stable for 24 months. Non-
hazardous. No MSDS required.

Note

For research use only

Isotype

Mouse / IgG1, kappa + Mouse / IgG1, kappa

Clonality

Monoclonal

MW

40-67kDa

Uniprot ID

[Q01546](#), [Q7Z794](#)

Entrez

[374454](#), [51350](#)

Dilution Range

Flow Cytometry (5ul per test per one
million cells or 5ul per 100ul of whole
blood);Immunofluorescence (1:50-
1:100);Optimal dilution for a specific
application should be determined.

Expiration Date

12 months from date of receipt.