

Product Datasheet

MLANA & SILV antibody (orb607976)

Description

Mouse monoclonal antibody to Melanoma

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant hMART-1 protein (DT101 & BC199); Extract of pigmented melanoma metastases from lymph nodes (HMB45)
Target	MLANA & SILV
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG2b, kappa + Mouse / IgG1, kappa
Clonality	Monoclonal
MW	20-22kDa (doublet) (MART-1): 90-100kDa (gp100)
Uniprot ID	P40967 , Q16655
Entrez	2315 , 6490
Dilution Range	Immunohistochemistry (Formalin-fixed) (1-2ug/ml for 30 minutes at RT)(Staining of formalin-fixed tissues requires heating tissue sections in 10mM Tris with 1mM EDTA, pH 9.0, for 45 min at 95°C followed by cooling at RT for 20 minutes); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.