

Product Datasheet

HPV-16 E6 + HPV-18 E6 antibody (orb607997)

Description

Mouse monoclonal antibody to HPV-16

	16/18)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	HPV18 E6- β -galactosidase fusion protein
Target	Not Applicable
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Mouse / IgG1, kappa
Clonality	Monoclonal
MW	16/17kDa
Uniprot ID	Not Applicable
Entrez	Not Applicable
Dilution Range	ELISA (For coating, order antibody without BSA); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.