

Product Datasheet

CD134 antibody (orb608014)

Description

Rat monoclonal antibody to CD134, Mouse

Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Purified recombinant Mouse OX-40 antigen
Preservatives	200ug/ml of Ab Purified from Bioreactor Concentrate by Protein A/G. Prepared in 10mM PBS with 0.05% rAlbumin & 0.05% azide. Also available WITHOUT rAlbumin & azide at 1.0mg/ml.
Concentration	Purified Ab with BSA and Azide at 200ug/ml
Storage	Antibody with azide - store at 2 to 8°C. Antibody without azide - store at -20 to -80°C. Antibody is stable for 24 months. Non-hazardous. No MSDS required.
Note	For research use only
Isotype	Rat / IgG1, kappa
Clonality	Monoclonal
MW	50kDa
Uniprot ID	P47741
Entrez	22163
Dilution Range	Does not block binding of OX40L; Flow Cytometry (0.5-1ug/million cells in 0.1ml); Immunofluorescence (1-2ug/ml); Optimal dilution for a specific application should be determined.
Expiration Date	12 months from date of receipt.