



Product Datasheet

Bax antibody (orb608054)

Description	Bax antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Bax
Target	Bax
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	21 kDa
Uniprot ID	Q07812
Dilution Range	WB=1:500-2000, ICC=1:100, IHC-P=1:100-500, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.



ICC analysis of SH-SY5Y cell using Bax m...



Western blot:Testis (Mouse);Testis (Rat)...