



Product Datasheet

PSMA1 antibody (orb608059)


Description	PSMA1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant human PSMA1
Target	PSMA1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	30 kDa
Uniprot ID	P25786
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human liver
carcinoma ...



IHC-P
analysis of
rat testis
using
PSMA1...



IHC-P
analysis of
mouse
cerebellum
using...