

Product Datasheet

GEF H1 antibody (orb608065)

Description	GEF H1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GEF H1
Target	GEF H1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	112 kDa
Uniprot ID	Q92974
Dilution Range	WB=1:500, IHC-P=1:50-200, ICC=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.