


Product Datasheet

NFKB P65 antibody (orb608066)

Description	NFKB P65 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human NFKB p65
Target	NFKB p65
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	61 kDa
Uniprot ID	Q04206
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:MCF-7(Human) Cell;Molt-4(Hu...