

Product Datasheet

SQSTM1 antibody (orb608068)

Description	SQSTM1 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human SQSTM1/P62
Target	SQSTM1/P62
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	60 kDa
CAS Number	Sequestosome-1
Uniprot ID	Q13501
Expiration Date	12 months from date of receipt.