



Product Datasheet

KLF4 antibody (orb608075)


Description	KLF4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Recombinant protein within human KLF4 aa 1-500
Target	KLF4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	56 kDa
Uniprot ID	O43474
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:100, Flow-Cyt=1ug/Test, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Flow cytometry of K562. Primary Antibody ...



ICC analysis of HeLa cell using KLF4 mon...



IHC-P analysis of Human esophageal cancer...