



Product Datasheet

HLA-DQA1 antibody (orb608079)


Description	HLA-DQA1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human HLA-DQA1
Target	HLA-DQA1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	26 kDa
Uniprot ID	P01906
Dilution Range	IHC-P=1:50-200, WB=1:500-2000, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
mouse colon
using HLA-...



IHC-P
analysis of
rat lymphoid
using HLA...



IHC-P
analysis of
human lung
carcinoma
u...