

Product Datasheet

PCNA antibody (orb608173)

Description	PCNA antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PCNA
Target	PCNA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	29 kDa
Uniprot ID	P12004
Dilution Range	IHC-P=1:200-800, IHC-F=1:200-800
Expiration Date	12 months from date of receipt.