



Product Datasheet

Beta Catenin antibody (orb608180)


Description	Beta Catenin antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human β -Catenin
Target	Beta catenin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	92 kDa
Uniprot ID	P35222
Dilution Range	IHC-P=1:200-1000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Human
descend
colon ca...



IHC-P
analysis of
rat colon
using Beta
C...



IHC-P
analysis of
mouse
kidney using
Bet...