

Product Datasheet

CD161 antibody (orb608186)

Description	CD161 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KLRB1. Please contact us for the exact immunogen sequence. The peptide is available as orb608188.
Target	CD161
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q12918
NCBI	002249, 002258
Entrez	3820
Dilution Range	IHC-P: 1:50~500, WB: 1:100~1000
Expiration Date	12 months from date of receipt.