

## Product Datasheet

### Tocilizumab (ACTEMRA) ELISA Kit (orb608202)

|                            |                                    |
|----------------------------|------------------------------------|
| <b>Description</b>         | Tocilizumab (ACTEMRA)<br>ELISA Kit |
| <b>Reactivity</b>          | Human                              |
| <b>Range</b>               | 0 ng/ml - 640 ng/ml                |
| <b>Tested Applications</b> | ELISA                              |
| <b>Note</b>                | For research use only              |
| <b>Application notes</b>   | Sandwich Assay                     |
| <b>Sample Types</b>        | Serum and Plasma                   |
| <b>Sensitivity</b>         | 10 ng/ml                           |
| <b>Expiration Date</b>     | Please enquire.                    |