



Product Datasheet

Piwil2 antibody (orb6086)

Description	Piwil2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Piwil2/Mili (885-973/973 aa)
Target	Piwil2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	107 kDa
Uniprot ID	Q8TC59
Dilution Range	WB=1:500-2000, Flow-Cyt=0.2µg /Test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot:testis (Mouse);293T Cell(Hu...)



Western blot:Testis (Mouse);Testis (Rat)...