

Product Datasheet

Phosphatase Dilution Buffer I [Out of stock] (orb608890)

Description	Phosphatase Dilution Buffer I [Out of stock]
Tested Applications	Phosphatase Assay
Target	Phosphatase Dilution Buffer I
Preservatives	Frozen solution of buffer containing: 25 mM HEPES pH 7.2, 50 mM NaCl, 2.5 mM EDTA, 5 mM DTT, 50 ng/ul rAlbumin.
Note	For research use only
Source	Stock solution of buffer used for diluting phosphatase protein.
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.