

Product Datasheet

GMCSF antibody (orb6090)

Description

GMCSF antibody

Species/Host

Rabbit

Reactivity

Guinea pig, Human, Porcine, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human GMCSF. Please contact us for the exact immunogen sequence. The peptide is available as orb1157.

Target

GM-CSF

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

16 kDa

Uniprot ID

[P48750](#), [P04141](#)

NCBI

[000749](#)

Entrez


[1437](#)

Dilution Range


WB: 1:100-1000, IHC-P: 1:100-600, IF/ICC: 1:100-600

Expiration Date


12 months from date of receipt.



IHC-P image of rat lung tissue using GMCSF...



IHC-P staining of pig small intestines tissue



Immunohistochemical staining of human tissue