


## Product Datasheet

### SLC2A4 antibody (orb609131)

|                            |                                                                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Goat polyclonal antibody to SLC2A4                                                                                                                          |
| <b>Species/Host</b>        | Goat                                                                                                                                                        |
| <b>Reactivity</b>          | Bovine, Canine, Human, Mouse, Porcine, Rabbit, Rat                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                |
| <b>Tested Applications</b> | ELISA, IF                                                                                                                                                   |
| <b>Target</b>              | SLC2A4                                                                                                                                                      |
| <b>Preservatives</b>       | Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing. |
| <b>Storage</b>             | Aliquot and store at -20°C. Minimize freezing and thawing.                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                       |
| <b>Application notes</b>   | Notes: Immunofluorescence: Strong expression of the protein seen in the membrane and cytoplasm of NIH3T3 cells. Recommended concentration: 10µg/ml.         |
| <b>Protein Sequence</b>    | NLEGPARKSLKRLT                                                                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                  |
| <b>MW</b>                  | 54.8                                                                                                                                                        |
| <b>Entrez</b>              | <b>6517</b>                                                                                                                                                 |
| <b>Dilution Range</b>      | P-ELISA: 1:16000, IF/ICC: 10 µg/ml                                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                             |



Immunofluorescence  
analysis of NIH3T3  
ce...