


Product Datasheet

GRIN2A antibody (orb609132)

Description	Goat polyclonal antibody to GRIN2A
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	GRIN2A (aa211-224)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the plasma membrane of NIH3T3 cells. Recommended concentration: 10µg/ml.
Protein Sequence	EDAKTQVQLKKIHS
Clonality	Polyclonal
MW	165; 144
Entrez	2903
Dilution Range	P-ELISA: 1:16000, IF/ICC: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH3T3 cells.