


Product Datasheet

MTOR antibody (orb609133)

Description	Goat polyclonal antibody to MTOR
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	MTOR / FRAP2 (aa1137-1150)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the ER and vesicles/cytoplasm of HEK293 cells. Recommended concentration: 10µg/ml.
Protein Sequence	ETVDRLTESLDFTD
Clonality	Polyclonal
MW	289
Entrez	2475
Dilution Range	P-ELISA: 1:128000, IF/ICC: 10 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH3T3 cells...