



Product Datasheet

Glucose 6 phosphatase alpha antibody (orb6097)

Description	Glucose 6 phosphatase alpha antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Glucose 6 phosphatase alpha (81-180/357 aa)
Target	G6PC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P35575
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=0.2ug/test
Expiration Date	12 months from date of receipt.



Paraffin-embedded human colon carcinoma ...



Paraffin-embedded rat kidney tissue fixe...