



Product Datasheet

GMCSFR beta antibody (orb6100)

Description	GMCSFR beta antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human GM-CSFR Beta (441-540/897 aa)
Target	IL3RB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	P32927
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human intestine cancer using GM...