



Product Datasheet

MVP Antibody (monoclonal, 8B12) (orb610911)


Description	MVP Antibody (monoclonal, 8B12)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human MVP recombinant protein (Position: A2-H259).
Form/Appearance	Lyophilized
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	100-110 kDa
Uniprot ID	Q14764
Dilution Range	Western blot, 0.1-0.5µg/ml, Human, Mouse, Rat Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat, By Heat Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10 ⁶ cells, Human
Expiration Date	12 months from date of receipt.



IF analysis of MVP using anti-MVP antibody.



Flow Cytometry analysis of A431 cells using MVP 8B12.



Flow Cytometry analysis of A549 cells using MVP 8B12.