



Product Datasheet

CD79a Antibody (monoclonal, 4G4) (orb610913)

Description	CD79a Antibody (monoclonal, 4G4)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	E.coli-derived human CD79a recombinant protein (Position: T121-P226). Human CD79a shares 91% amino acid (aa) sequence identity with mouse CD79a.
Form/Appearance	Lyophilized
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	44 kDa
Uniprot ID	P11912
Dilution Range	Western blot, 0.1-0.5µg/ml, Human Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat
Expiration Date	12 months from date of receipt.



IHC analysis of CD79A using anti-CD79A
a...



WB analysis using anti ZO-1 antibody.Lan...