


Product Datasheet

Glucocorticoid Receptor (phospho-Ser211) antibody (orb6110)


Description

Glucocorticoid Receptor (phospho-Ser211)

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Glucocorticoid Receptor around the phosphorylation site of Ser211 (NE(p-S)PW)
Target	Phospho-Glucocorticoid Receptor (Ser211)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85 kDa
Uniprot ID	P04150
Dilution Range	Flow-Cyt=1ug/Test, WB=1:500-2000, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of A549.Primary Antibody ...



Flow cytometry of Mouse spleen (Black).P...