

Product Datasheet

Human TRIM21 Protein (Biotin) [Out of stock] (orb611144)

Description	TRIM21 is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding...
Species/Host	Insect
Immunogen	TRIM21 is expressed with a -6x His tag at N-terminus
Preservatives	TRIM21 solution is supplied in 20mM HEPES pH-7.6, 0.01mM EDTA and 0.02% SDS
Form/Appearance	Sterile Filtered clear solution
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Binds IgG type human auto antibodies. Functional Streptavidin based ELISA test (analysis of positive/negative samples.)
Purity	Greater than 80% as determined by SDS-PAGE
MW	TRIM21 Human Recombinant, Biotin produced in SF9 is a glycosylated, polypeptide chain having a calculated molecular mass of 52kDa
Expiration Date	6 months from date of receipt.