

Product Datasheet

Bovine TROVE2 Protein (orb611146)

Description	Recombinant Bovine TROVE2 protein
Preservatives	TROVE2 is supplied in 20mM HEPES buffer pH-8.0, 150mM NaCl and 20% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Purity	Greater than 83.0% as determined by SDS-PAGE.
Source	Bovine tissues
Expiration Date	12 months from date of receipt.