

Product Datasheet

5-Propargylamino-dCTP-XX-Cy3 (orb611225)

Description	5-Propargylamino-dCTP-XX-Cy3
Form/Appearance	solution in 10 mM Tris-HCl; Color: pink; pH: 7.5 ± 0.5
Concentration	1.0 mM-1.1 mM
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 550 nm, λ_{em} 570 nm, ϵ 150.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
Formula	C ₅₅ H ₇₇ N ₈ O ₂₂ P ₃ S ₂
Purity	≥ 95% (HPLC)
MW	Theoretical MW: 1359.30 g/mol (free acid); Detected MW: 1358.38 g/mol (free acid)
SMILES	<chem>N(CC#Cc1cn(c(=O)nc1N)[C@@H]1O[C@H](COP(=O)(O)OP(=O)(O)OP(=O)(O)O)[C@H](C1)OC(=O)CCCCCNC(=O)CCCCCNC(=O)CCCCCN1c2ccc(S(=O)(=O)[O-])cc2C(C)(C)/C1=C\C=C1=*(c2ccc(cc2C1(C)C)S(=O)(=O)O)CC</chem>
Expiration Date	12 months from date of receipt.