

Product Datasheet

HighYield T7 Cy3 RNA Labeling Kit (orb611239)

Description	HighYield T7 Cy3 RNA Labeling Kit (UTP-based) is designed to produce randomly Cy3-modified RNA probe...
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 550 nm, λ_{em} 570 nm, ϵ 150.0 L mmol ⁻¹ , cm ⁻¹ (Tris-HCl pH 7.5).
SMILES	<chem>P(=O)(OP(=O)(O)O)(O)OP(=O)(O)OC[C@H]1O[C@H]([C@@H]([C@@H]1O)O)n1c(=O)[nH]c(=O)c(c1)/C=C/CNC(=O)CCCCCNC(=O)CCCCCN1c2cc c(S(=O)(=O)O)cc2C(C)(C)/C/1=C\C=C\C1=[N+](c2ccc(cc2C1(C)C)S(=O)(=O)[O-])CC</chem>
Hazard Information	Non-Toxic
Expiration Date	12 months from date of receipt.