

Product Datasheet

Human Siglec-9 Protein (orb612118)

Description

Human Siglec-9, Fc Tag (SI9-H5256) is

It contains AA Gln 18...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Siglec-9

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 62.5 kDa. The protein migrates as 80-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_055256.1

Purity

● -20°C to -70°C for 12 months in lyophilized state; ● -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Source

Human Siglec-9, Fc Tag (orb612118) is expressed from human 293 cells (HEK293). It contains AA Gln 18 - Gly 348 (Accession # Q9Y336-1).

MW

62.5 kDa

Uniprot ID

[Q9Y336](#)

NCBI

[NP_055256.1](#)

Expiration Date

6 months from date of receipt.