

Product Datasheet

Human Siglec-9 Protein (orb612120)

Description	siglec-9 (HGMW-approved symbol SIGLEC9) a member of the sialic acid-binding Ig-like lectin (Siglec) ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Siglec-9
Preservatives	25 mM MES, 150 mM NaCl, pH5.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 37.9 kDa. The protein migrates as 53-66 kDa under reducing (R) condition (SDS- PAGE) due to glycosylation.
Protein Sequence	NP_055256.1
Purity	90%
Source	Human Siglec-9, His Tag (orb612120) is expressed from human 293 cells (HEK293). It contains AA Gln 18 - Gly 348 (Accession # Q9Y336-1).
MW	37.9 kDa
NCBI	NP_055256.1
Expiration Date	6 months from date of receipt.