

Product Datasheet

Human IL-13 R alpha 2 Protein (orb612134)

Description Interleukin-13 receptor subunit alpha-2 is also known as IL13R α 2, IL13Ra2 cluster of differentiatio...

Reactivity Human

Endotoxins 1.0 EU per μ g

Conjugation Unconjugated

Target IL-13 R alpha 2

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag C-10 $\times$ His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 39.0 kDa. The protein migrates as 45-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.


Purity 90%

Source Human IL-13 R alpha 2, His Tag (orb612134) is expressed from human 293 cells (HEK293). It contains AA Asp 27 - Arg 343 (Accession # Q14627-1).


MW 39.0 kDa

NCBI [NP_000631.1](#)


Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human IL-13 R alpha...



Line graph illustrates about the Ag-Ab r...



Line graph illustrates about the Ag-Ab r...