



Product Datasheet

GRP78 antibody (orb612173)

Description	GRP78 antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human GRP78/Bip
Target	GRP78/Bip
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	78 kDa
Uniprot ID	P11021
Dilution Range	WB=1:500-1000, ICC=1:100
Expiration Date	12 months from date of receipt.



Western blot:U87MG(Human) Cell;Hela(Huma...



ICC analysis of HepG2 cell using GRP78/B...