


Product Datasheet

Hydrogen Exchanger 3 antibody (orb612180)

Description	Hydrogen Exchanger 3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human hydrogen exchanger 3 (581-680/834 aa)
Target	hydrogen exchanger 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93 kDa
Uniprot ID	P48764
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Kidney (Mouse); Kidney (Rat)...