


Product Datasheet

MNX1 antibody (orb612185)

Description	MNX1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human MNX1/HLXB9 (131-230/401 aa)
Target	MNX1/HLXB9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P50219
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Cerebrum(Mouse);Cerebrum(Ra...)