


Product Datasheet

GASC1 antibody (orb612193)

Description	GASC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human GASC1/KDM4C (291-390/1056 aa)
Target	GASC1/KDM4C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120 kDa
Uniprot ID	Q9H3R0
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:A431 (Human) Cell;HepG2 (Hu...)