

Product Datasheet

Chondroitin Sulfate antibody (orb612200)

Description	Chondroitin Sulfate antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Chondroitin sulfate sodium salt from bovine cartilage (purified)
Target	Chondroitin Sulfate
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
CAS Number	39455-18-0
Expiration Date	12 months from date of receipt.