



Product Datasheet

PCNA antibody (orb612208)


Description	PCNA antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant human PCNA protein
Target	PCNA (Nuclear Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	29 kDa
Uniprot ID	P12004
Dilution Range	WB=1:1000-3000, IHC-P=1:100-500, Flow-Cyt=1ug/Test, IHC-F=1:100-500, ICC=1:50, IF=1:100-500
Expiration Date	12 months from date of receipt.



Flow cytometry of HeLa. Primary Antibody ...



IHC-P analysis of mouse testis using PCN...



IHC-P analysis of rat stomach using PCNA...