


Product Datasheet

TNFR2 antibody (orb612211)

Description	TNFR2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human TNFR2
Target	TNFR2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	46 kDa
Uniprot ID	P20333
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: MCF-7 Cell; Lane 2: Ju...