

Product Datasheet

WNT4 antibody (orb612217)

Description	WNT4 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human WNT4 protein (30-350/350 aa)
Target	WNT4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37 kDa
Uniprot ID	P56705
Expiration Date	12 months from date of receipt.