


Product Datasheet

Galectin 8 antibody (orb612219)

Description	Galectin 8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Galectin 8
Target	Galectin 8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	35 kDa
Uniprot ID	O00214
Dilution Range	IHC-P=1:50-200, WB=1:500-1000, IHC-F=1:50-200, ICC=1:400-800, IF=1:400-800
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human liver
carcinoma ...