



Product Datasheet

VEGFR1 antibody (orb612220)

Description	VEGFR1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human VEGFR1
Target	VEGFR1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	147 kDa
Uniprot ID	P17948
Dilution Range	WB=1:1000-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of human brain using VEGF...



Western blot: Lane 1: Placenta (Mouse); Lan...