



Product Datasheet

CD71 antibody (orb612221)


Description	CD71 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Transferrin receptor
Target	TFRC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	95 kDa
Uniprot ID	Q62351
Dilution Range	WB=1:500-2000, IHC-P=1:50-200, ICC=1:100, IHC-F=1:50-200
Expiration Date	12 months from date of receipt.



ICC analysis of HeLa cell using TFRC mon...



IHC-P analysis of rat placenta using CD7...



IHC-P analysis of human lung carcinoma u...