



Product Datasheet

Tubulin Beta-III antibody (orb612226)

Description	Tubulin Beta-III antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Tubulin beta-III
Target	Tubulin beta-III
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	50 kDa
Uniprot ID	Q13509
Dilution Range	WB=1:1000-2000, IHC-P=1:50-200, ICC=1:50, Flow-Cyt=1:50, IHC-F=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: PC-12 Cell; Lane 2: HeLa



Western blot: Lane 1: Mouse Cerebrum; Lane 2: ...