


Product Datasheet

CD4 antibody (orb612230)

Description	CD4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human CD4
Target	CD4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	48 kDa
Uniprot ID	P01730
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/Test, IHC-P=1:100-500, IHC-F=1:50-200, ICC=1:50-200, IF=1:50-200
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: U937 (Human) Cell; La...